










# TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



## PRODUCT ADVANTAGES



EABS Brake system

500W dual motor

High battery life

360° joystick



Hand/electric mode

Easy to fold

Aluminum alloy frame

Anti-backover design

## PRODUCT PARAMETER

| Parameter                  | Category                         |
|----------------------------|----------------------------------|
| Series/Model No.           | Alphard Series/74501             |
| Frame material             | Aluminum alloy                   |
| Unfolded Dimension         | 1000*610*930mm                   |
| Folded Dimension           | 660*380*770mmmm                  |
| Motor                      | DC 24V、250W                      |
| Battery                    | Lithium battery DC 24V 13 AH     |
| Weight                     | 24KG                             |
| Loading Capacity           | 120KG                            |
| Front wheel specifications | 8"*2" (200*50mm) PU tire         |
| Rear wheel specifications  | 12"*2" (305*60mm) pneumatic tire |
| Max Speed                  | 6KM/H                            |
| Driving Range              | 15KM                             |

## PRODUCT STRUCTURE



Motor 1

Motor 2

500W dual motor,  
DC motor /strong climbing ability



LG cell lithium batteries  
travel distance about 15km

## PRODUCT DETAILS



### FOLDING DESIGN

Lightweight folding portable design  
Can be placed in the trunk, Similar to suitcase size.



### LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



### THE ARMREST CAN BE RAISED

Handrails can swing up and down  
Easy to hold the patient up and down wheelchair



### HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel

