

Hersteller von elektrischen Rollstühlen in China

TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



PRODUCT ADVANTAGES

| | | | |
|--------------------|-----------------|----------------------|----------------------|
| | | | |
| EABS Brake system | 500W dual motor | High battery life | 360° joystick |
| | | | |
| Hand/electric mode | Easy to fold | Aluminum alloy frame | Anti-backover design |

PRODUCT PARAMETER

| Parameter | Category |
|----------------------------|----------------------------------|
| Series/Model No. | Alphard Series/74501 |
| Frame material | Aluminum alloy |
| Unfolded Dimension | 1000*610*930mm |
| Folded Dimension | 660*380*770mm |
| Motor | DC 24V、250W |
| Battery | Lithium battery DC 24V 13 AH |
| Weight | 24KG |
| Loading Capacity | 120KG |
| Front wheel specifications | 8"*2" (200*50mm) PU tire |
| Rear wheel specifications | 12"*2" (305*60mm) pneumatic tire |
| Max Speed | 6KM/H |
| Driving Range | 15KM |

PRODUCT STRUCTURE



Motor 1
500W dual motor,
DC motor /strong climbing ability



Motor 2
LG cell lithium batteries
travel distance about 15km

PRODUCT DETAILS



FOLDING DESIGN

Lightweight folding portable design
Can be placed in the trunk, Similar to suitcase size.



LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



THE ARMREST CAN BE RAISED

Handrails can swing up and down
Easy to hold the patient up and down wheelchair



HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel

