

Produsen kursi roda listrik Cina

TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



PRODUCT ADVANTAGES

EABS Brake system	500W dual motor	High battery life	360° joystick
Hand/electric mode	Easy to fold	Aluminum alloy frame	Anti-backover design


PRODUCT PARAMETER

Parameter	Category
Series/Model No.	Alphard Series/74501
Frame material	Aluminum alloy
Unfolded Dimension	1000*610*930mm
Folded Dimension	660*380*770mm
Motor	DC 24V、250W
Battery	Lithium battery DC 24V 13 AH
Weight	24KG
Loading Capacity	120KG
Front wheel specifications	8"*2" (200*50mm) PU tire
Rear wheel specifications	12"*2" (305*60mm) pneumatic tire
Max Speed	6KM/H
Driving Range	15KM

PRODUCT STRUCTURE



Motor 1
500W dual motor,
DC motor /strong climbing ability



Motor 2
LG cell lithium batteries
travel distance about 15km

PRODUCT DETAILS



FOLDING DESIGN

Lightweight folding portable design
Can be placed in the trunk, Similar to suitcase size.



LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



THE ARMREST CAN BE RAISED

Handrails can swing up and down
Easy to hold the patient up and down wheelchair



HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel

