




Lifting from floor



Toileting/Bathing transfer



Bed-Chair transfer

Patient lift is an assistive device, suitable for transfer people such as Wheelchair Users. Bariatric Patient. Bedridden patients, Disabled, Handicap, lower limb fracture patients, elderly people with poor ability to walk to a wheelchair.

This electric battery power patient lift offers exceptional value in an electric patient lifter. It's an affordable powered lift with battery and charger.

Lift from floor, bed, commode, or wheelchair. Easy to use Lift system has a wide cradle that accommodates a 2-point or 4-point sling.

## Features:

- Heavy-gauge construction, with 6-Point Cradle
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Fast, easy charge from AC outlet
- To charge, simply plug in unit. Battery does not need to be removed to be charged
- Lift will not operate while plugged in
- For home-care use
- Optional Weight Scale

|                        |   |
|------------------------|---|
| <b>Product Name</b>    | Battery Powered Patient Lift                                |
| <b>Net Weight</b>      | 49kg  |
| <b>Weight Capacity</b> | 204kg   |
| <b>Height</b>          | 27.5"-76"   |
| <b>Width of Base</b>   | 24"-41"   |
| <b>Function</b>        | Handicapped or disabled transfer, walking assist, home care |
| <b>Application</b>     | Home care, nursing home or hospital, etc.                   |
| <b>Thrust</b>          | 10000N  |

## Detailed Description:



**Easy operating**

During the lifting process, the hanger can be rotated 360 to obtain a comfortable position. Optional weight scale.



**Emergency stop button**

with Emergency stop button and emergency lowering system. The pump can be switched to manual mode to lower the patient safely.



**Grip handles**

Large, easy to grip, for easy maneuverability.



**Castors with brake**

Brakes in the rear wheels prevent the device from moving during use.