

This is a quiet and hard-working compressor nebulizer. It offers higher flow flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.  
**Detailed Description:**  
Ultra compact design makes it extremely portable.  
Powerful piston pump is efficient, dependable and quiet.  
Delivers consistent particle size of 5.5 microns and fast treatment times.



Suitable for A Variety of Situations



**PORTABLE COMPRESSOR NEBULIZER**


Penetrate into Bronchus & Lung  $2\text{--}5\mu\text{m}$



Adults & Children

**MICRON NEBULIZATION FINE PARTICLES**

MMAD:  $3.6 - 4.0 \mu\text{m}$




**QUIET ENVIRONMENT SILENT ATOMIZATION**


Low noise design, easy for children to sleep



**USE BASIC STEPS**




More Details



**Micron level nebulisation—2µm—5µm**  
Nebulised particles directly reach the upper respiratory tract, trachea, alveoli and lower respiratory tract, making it easier to absorb the medication.



**High atomisation efficiency**  
Smooth atomisation for continuous operation.



**Low noise design**  
Low power, low vibration, baby atomisation does not resist.