

Clamp-On Tub Rail

The height adjustable bathtub grab bar safety bar is suitable for almost any style of bathtub. The attractive white powder coated steel frame is both durable and beautiful. The grab bar is adjustable from 14" to 17", allowing you to find the perfect height for your needs. The grab bar is angled parallel to the shower wall. Adjustable width clips are attached by a swivel locking mechanism that locks the rubber mat to all tub contact points. The rubber pad prevents scratches to the tub surface and ensures a secure fit. The bathtub grab bar safety bar is easy to install and requires no tools. Please note that this safety bar is not suitable for use with fiberglass bathtubs.



Product Advantage

- 14.5" - 17" height adjustment (for RTL15020 only)
- Width adjusts from 3" - 7"
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch-proof security
- Fits most tubs (not for use with fiberglass tubs)



Durable, white powder-coated steel construction, Tool-free installation;
Securely tightens to fit edges 3" - 7" in width