

Three Motors Full Electric Low Bed

Quiet, convenient, full electric low height bed provides smooth position change with easy-to-use Hand Pendant control, and overall height of the bed frame for maximum comfort and versatility for the patient or caregiver.



Electric Beds Buying Guide



Weight Capacity 204kg



Low Bed



Castors with Brake



Full Electric Bedc
with Three Motors



Product Structure

- ① Foot Spring
- ② Head Spring
- ③ Head/Foot Board
- ④ Emergency Manual Hand Crank
- ⑤ Casters
- ⑥ Hand Pendant Control
- ⑦ 3 Motors/Actuators

Product Advantage

- 1.The foot section positioning and overall height of the bed frame for maximum comfort and versatility for the patient or caregiver
- 2.Channel frame construction provides superior strength and reduced weight
- 3.Reinforced frame resists twisting and bending
- 4.Long lasting Zinc-coated spring deck for Patient comfort
- 5.Heavy-duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- 6.Headboard is taller than footboard to blend with home décor
- 7.Labels and color-coded springs help ensure accurate installation of side rails for patient safety
- 8.High performance workhorse muscle of the 3-motor power plant is an industry standard
- 9.In the event of power failure, top mounted Hi/Lo crank accessibility and manual crank motor override of sleep surface adjustments