

Manual Patient Mover 72110

SPECIFICATIONS

Size (L*W*H)	35.8" * 25.2" * 4.2"
Knee Pad Measurements	16.5" * 9.8"
Handle Height (from base)	35.4"
Weight	74lbs
Maximum User Weight Limit	441lbs
Castors	4"
Ground Clearance	5.3"
Underbed Clearance	5.9"




441 lbs
weight capacity

DETAILED DESCRIPTION




Flip up the Seat Pads so the User can occupy the Able Assist□

Support provided by knee pads for users;





You can use the able assist in the bathroom and can be wheeled up to a toilet , the wheels can be positioned at either side of the toilet offering the user more dignity when sitting down.





A substantial aluminium robust frame, The Able Assist Transfer Aid can be used around the home to help provide a safer and easier way of transferring a client. Maximum user weight of 200kg.



Fitted with 4 wheels that enhance manoeuvrability and enable users to be transferred with ease, 2 of the wheels have breaks that ensure maximum security and stability when transferring a patient.

SIMILAR PRODUCTS

