

Procyon compressor nebulizer provides fast, efficient and reliable nebulized inhalation therapy, designed for home use. This is a quiet and hard-working compressor nebulizer suitable for nebulizing most aerosol medications prescribed for common respiratory diseases. It offers higher liter flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

Detailed Description: Makes it extremely portable.

It is compact, portable, a efficient, dependable and quiet.

Powerful, gentle, quiet & efficient nebulizer and quiet.

Delivers consistent particle size of 4.5 microns and fast treatment times.

Series / Model No.	Procyon series / 70000
Unit Dimensions	191mm x 102mm x 117mm
Weight	4.1kg (11.04kg)
Electrical Requirements	100-240VAC, 50/60Hz
Operating Noise Level	< 60 dBA
Maximum Pressure	< 15psi
Compressor flow air flow	> 6L/min
Particle Size	2.0 - 5.0µm
Nebulizer Flow	> 2.0L/min
Residual Liquid	10.5ml




PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



MICRON NEBULIZATION FINE PARTICLES

MMAD: 3.6 - 4.0 µm




QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep




USE BASIC STEPS



Suitable for A Variety of Situations




More Details



Micron level administration — 2µm–5µm
Nebulized particles directly reach the upper respiratory tract, trachea, bronchi and lower respiratory tract, making it easier to absorb the medication.



High atomization efficiency
Special atomization for continuous operation.



Low noise design
Low power, low vibration, baby atomization does not react.