

This is a quiet and hard-working compressor nebulizer. It offers higher liter flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

**Detailed Description:**


Ultra compact design makes it extremely portable.

Powerful piston pump is efficient dependable and quiet.


Delivers consistent particle size of ≤5 microns and fast treatment times.



**Suitable for A Variety of Situations**



## More Details




### **Micron level nebulization—0.5 $\mu$ m~5 $\mu$ m**

Nebulized particles directly reach the upper respiratory tract, trachea, alveoli and lower respiratory tract, making it easier to absorb the medication.



### **High atomization efficiency**

Smooth atomization for continuous operation.



### **Low noise design**

Low power, low vibration, baby atomization does not resist.

