

# Power Hoyer Lift 71920



Weight Capacity  
450lbs. / 204kg



Emergency  
Stop Button



6-Point Cradle



Battery Control Box



Optional Weight Scale



## Product Description

This Electric Battery Power Patient Lift offers exceptional value in an Electric Patient Lifter. It's an affordable

powered lift with battery and charger.

Lift from floor, bed, commode, or wheelchair. Easy to use Lift system has a wide cradle that accommodates

a 2-point or 4-point sling.


White powder-coated finish is easy to clean and won't stain or tarnish.

## Features:

1. Heavy-gauge construction
2. Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
3. Fast, easy charge from AC outlet
4. To charge, simply plug in unit Battery does not need to be removed to be charged Lift will not operate while plugged in For home-care use
5. 6-point swivel attachment allows for 2, 4 or 6 point slings
6. Base height only 4.75" for low under bed clearance
7. Multiple sling options available
8. Lifts up to 450 lbs.
9. Optional Weight Scale.

## Applications:

## Detailed Description



Lifting from floor



Bed-Chair transfer



Toileting/Bathing

Lift from floor, bed, commode, or wheelchair.

Easy to use Lift system has a six point cradle that will accommodate a six point sling, 4 point sling or a 2 point slings.

Optional Weight Scale



With Emergency stop button and emergency lowering system.

This Lift has a removable battery pack that plugs into the control box/charger.

The pump even includes an emergency button and can be switched to manual mode to lower the patient safely.



Adjustable base,Expand base of lifter to widest possible position to maximize stability and ease of use.

## Product Comparison:



	Battery Powered Patient Lift 71920	Battery Powered Patient Lift 71921	Bariatric Battery Powered Lift 71922	Electric Patient Lift71940
 Weight Capacity	450 lbs/204kg	450 lbs/204kg	600 lbs/272kg	400 lbs/181kg
 Net Weight	49kg	47.6kg	59.9kg	39.9kg
 Gross weight	55.3kg	54.4kg	67.1kg	48.5kg
 Base Height	4.75"	4.5"	4.5"	4.4"
 Boom Height	27.5" - 76"	29" - 80.5"	32" - 75.5"	20" - 65"
 Casters	3"Front/5"Rear	3"Front/4.5"Rear	3" Front/4" Rear	4"
 Min and Max Base Width	41"/24"	41"/24"	45.75" / 25"	33" / 23"
 Weight Scale	Optional	/	Optional	/