

Power Patient Hoist 71920

SPECIFICATIONS

Boom Height	27.5"-76"
Base Width	41" (open) / 24" (closed)
Base Height	4.75"
Casters	3" Front/5" Rear
Weight	49kg
Carton Shipping Weight	55.3kg
Weight Limit	450 lbs/204kg
Product Material	Steel



450 lbs
weight capacity

APPLICATIONS

DETAILED DESCRIPTION



Lift from floor, bed, commode, or wheelchair.

Easy to use Lift system has a six point cradle that will accommodate a six point sling, 4 point sling

or a 2 point slings.

Lifting from floor



Bed-Chair transfer



Toileting/Bathing



with Emergency stop button and emergency lowering system.


This Lift has a removable battery pack that plugs into the control box/charger.

The pump even includes an emergency button and can be switched to manual mode to lower the patient safely.



Wide base provides additional stability,

Base height (under bed clearance) is 4.75" for at home patient lift use.





Optional Weight Scale.

SIMILAR PRODUCTS

