

Power Patient Hoist 71940

SPECIFICATIONS

Boom Height	52"
Weight Limit	400lbs/181kg
Base width	23"/33"
Base Height	4.4"(110mm)
Weight	39.9kg
Lifting range	20"-65"
Casters	4"
Voltage Rating	24V DC



400 lbs
weight capacity

APPLICATIONS



Lifting from floor



Bed-Chair transfer



Toileting/Bathing

DETAILED DESCRIPTION



Lift from floor, bed, commode, or wheelchair.


Easy to use Lift system has a six point cradle that will accommodate a six point sling, 4 point sling or a 2 point slings.




with Emergency stop button and emergency lowering system.


This Lift has a removable battery pack that plugs into the control box/charger.

The pump even includes an emergency button and can be switched to manual mode to lower the patient safely.





A substantial steel robust frame,Maximum user weight limit of 181Kg.It is very manoeuvrable and light to use





The strong sturdy base. Rear casters adjust easily and lock securely so the lift will not move when desired.

SIMILAR PRODUCTS

