

Sit to Stand Transfer Aid

The Sit to Stand Transfer Aid assists with pivoting transfers between different seat surfaces and can transport Users short distances.Reduces strain on carer and enables safe and comfortable turning transfers for people who require additional support when standing.



Sit to Stand Transfer Aid Buying Guide



Weight capacity 150kg




Adjustable leg support



Sturdy steel frame



with suspenders



Product Structure

- Handle
- Steel frame
- Knee pads
- Knob
- Foot pedal
- Wheels

Product Advantage

- Its narrow, recessed front fits between wheelchair wheels and around toilets for perfect positioning and padded lower leg supports can be adjusted in width, angle and height to suit user comfort.
- Foot operated brakes on centre wheels provide stability during standing/sitting movement.
- Substantial base area provides stability. Easily dismantled/reassembled without tools for transportation and storage.