

# Ultra Light Full Electric Low Bed

## SPECIFICATIONS

Bed Height Minimum Without Casters	9.5"
Bed Height Maximum With Casters	23.5"
Unit Dimensions	36" (W) x 88" (L)
Weight	175 lbs.
Carton Shipping Weight	191 lbs.
Weight Limit	450 lbs.



**450 lbs**  
weight capacity

## DETAILED DESCRIPTION



This type care bed is especially for homecare sector.  
Superior Strength,High-strength frame provides full 450-lb weight capacity.



The Full-Electric Low Hospital Bed provides safety and convenience for a Home Care.  
All positions are controlled electrically to provide greatest convenience for the patient and ease to the caregiver.



The hand pendant controls the positioning of the upper body and or knees, and can be placed anywhere for maximum patient comfort.  
The pendant control also adjusts bed-frame height for efficient use of accessories and for safe and easy transfer.



Low deck height, This Low Hospital Bed has 9.5" to 23.5" deck height. The split-spring design incorporated into the low bed profile is designed to prevent injuries from falling out of bed.

## SIMILAR PRODUCTS

