


Home Bed Assist Handle

For standard Home Style Beds, this hand rail has a support bar that goes between mattress and box springs, and an adjustable strap to secure it. Can be mounted on either side of bed and provides a stable support. The removable power-coated steel handle includes a mid bar, which creates a reliable grip at any height and can be detached when not in use.



Product Advantage

Assistance for getting in and out of bed. Can be used on either side of the bed.

- Removable handle.
- Easy, quick, no tool assembly. Comes in two pieces.
- Safety strap wraps around mattress or box spring to ensure a secure fit.
- Special extra-long safety strap enables user to be in a single, twin, double, queen or king size bed.
- Made of powder coated steel for long life and easy cleaning.
- Handle with multiple crossbars creates easy grip for any height use.



The home bed-assist rail provides patients with assistance getting into and out of the bed. This model features a safety strap you wrap around the mattress or box spring to ensure a safe and secure fit.