

Portable Commode Wheelchair

The Commode Wheel Chair is perfect for people who need a commode but also want the convenience of a wheelchair for easy transportation and mobility. The rear wheels are equipped with brake casters to prevent sliding after sitting down. The seat cut-out provides comfort and allows the chair to be positioned over a standard toilet for use, but also comes with a bucket to be used anywhere it is required.



Commode Wheelchair Buying Guide



Weight Capacity
(350 lbs)



Safety
Locked casters



Material



Convenience



Product Structure

- Back Support
- Armrest
- Cover/cushion
- Bucket
- Swivel casters (2 with locks)
- Footrests

Product Advantage

- Chrome-plated-coated steel tubing frame with welded construction□also available in silver vein finish.
- Comes standard with padded seat.
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Swing-away and detachable footrests .
- Injection-molded large plastic toilet seat, 12 quart bucket and cover come standard.
- Removable, tool-free back.
- Rust-resistant, swivel casters (2 with locks).
- Folding, Portable, Upholstered Commode.