

Ultra-Light Electric Bed

A great Ultra-Light Electric Bed for your personal care!
Ultra-Light Electric Bed is ideal for any individual who needs a bed at home, in a hospital, or at a nursing home. A self-contained motor remains incredibly quiet, and motor assembly can be installed and removed while patient remains in bed.
Comfort, Convenience, Quality : Easily adjust the bed to the height and position you need. The bed frame is simple to assemble, cleans easily.



Electric Bed Buying Guide



Weight Capacity
(450lbs)



Functions and features



Material



Accessories



Product Structure

- Guardrail (Optional)
- Motor
- Spring bed
- Head/Footboard
- Lock pin
- Casters


Product Advantage

The bed features channel frame construction with zinc coated spring deck to provide superior strength and reduced weight. The bed is strong enough to withstand a weight capacity of up to 450 pounds yet features a caster that moves easily, making adjustments between 9.5"-20" deck height a cinch.

Comes with 9V battery to lower the bed section up to 9 times during power outage.

Headboard is taller than footboard to better fit home decor.

Bed accessories



Hand Control Pendant

6-button control with glow in the dark buttons.

Even with the lights out the pendant can be quickly located by the user or the attendant.

Compatible with all bed mattresses

Any standard size mattress will fit this frame. Whether you want to use a foam, spring or alternating pressure air mattress, this frame will be perfect.

Half Rail Included:

These safety rails clamp on and off easily. An easy release button enables each rail to be independently lowered or raised simply.