

# Ultra-Light Electric Bed

A great Ultra-Light Electric Bed for your personal care!  
Ultra-Light Electric Bed is ideal for any individual who needs a bed at home, in a hospital, or at a nursing home. A self-contained motor remains incredibly quiet, and motor assembly can be installed and removed while patient remains in bed.  
Comfort, Convenience, Quality : Easily adjust the bed to the height and position you need. The bed frame is simple to assemble, cleans easily.



## Electric Bed Buying Guide



Weight Capacity  
(450lbs)



Functions and features



Material



Accessories



## Product Structure

- Guardrail (Optional)
- Motor
- Spring bed
- Head/Footboard
- Lock pin
- Casters


## Product Advantage

The bed features channel frame construction with zinc coated spring deck to provide superior strength and reduced weight. The bed is strong enough to withstand a weight capacity of up to 450 pounds yet features a caster that moves easily, making adjustments between 9.5"-20" deck height a cinch.

Comes with 9V battery to lower the bed section up to 9 times during power outage.

Headboard is taller than footboard to better fit home decor.

Bed accessories



**Hand Control Pendant**

6-button control

Even with the lights out the pendant can be quickly located by the user or the attendant.

**Compatible with all bed mattresses**

Any standard size mattress will fit this frame. Whether you want to use a foam, spring or alternating pressure air mattress, this frame will be perfect.

**Half Rail Included:**

These safety rails clamp on and off easily. An easy release button enables each rail to be independently lowered or raised simply.