

HOME BED ASSIST RAIL



The home bed-assist rail and bed board combo provides patients with assistance getting into and out of the bed. This model features a tri-fold bed board for easy portability and storage that simply slides under the mattress for extra firm support. The removable power-coated steel handle includes a "mid bar," which creates a reliable grip at any height, and can be detached when not in use. The product arrives in two pieces and requires simply, quick tool-free assembly.

PRODUCT ADVANTAGE

- Assistance for getting in and out of home style bed
- Can be used on either side of the bed
- Tri fold bed board for easy portability and storage
- Removable handle
- Removable powder coated steel handle enables user to remove when not in use
- Easy, quick, no tool assembly. Comes in 2 pieces
- Simply slides under the mattress for extra firm support
- Handle with the "Mid Bar", creates grip for any height use

The 75040 is ideal for limited mobility

individuals who require assistance getting into and out of bed.
Tool-free assembly - Can be assembled without difficulty in the comfort of your own home, no tools required
Removable handle - Store the handle when not in use for a virtually undetectable look
Extra firm design - Bed board uses weight of mattress to remain secure while handle is in use, meaning the unit remains securely in place with no permanent installation required
Folding Bed Board - Conveniently folds into three parts for easy storage and transportation

