

## Hydraulic patient lift

Deneb series lifts are very easy to maneuver and offer a wide selection of slings for a variety of lifting situations. Transfers can be to and from beds, chairs, floor to bed, lateral transfers, bathroom and toilet. The equipment has a compact and versatile design, suitable for moving all types of patients. Light and easy to disassemble, convenient for quick transportation and installation. It can make patients and caregivers feel safe, comfortable and stable.



### Patient Lift Buying Guide



weight capacity  
(450 pounds)




Lift Size Choose



Six point rotating bar



Hydraulic pump



**product structure**

- Hanger
- slingshot
- Basic
- rear wheel
- Base adjustment lever
- Hydraulic pump
- Pen
- push hand

**product advantage**

**Patient Elevator (Hoyer Elevator) Features:**

- The lift features a range of 20"-64" and allows people to be lifted from a lying position on the floor
- Pump Handle: Able to swivel from side to side for caregiver comfort. How many times can it be used?
- Tips: It is highly recommended that patients be assisted by two people while lifting.
- Easy assembly, it only takes you 2 minutes to assemble the lifts.

**Lift Size Choose**



**Sling Features:**

Easy to Use - The sling helps transfer users from a bed to a wheelchair, recliner, or shower chair. The large opening in the design also allows you to go to the bathroom.

Maintenance (easy to clean): Washable polyester fabric can be used in wet or dry environments.

Multiple Use: Ideal for transferring a patient from bed to chair, wheelchair or commode, chair to chair or from floor to bed or chair using your manual or electric lift.