

Rollator Walker

The rollator restores activity and independence for the elderly and mobility challenged. Is a light and versatile rollator which is approved for indoor and outdoor use. With a width of only 27", the Compact rollator is even easier to manoeuvre in tight spaces, allowing you to go where you want when out on the town or when travelling. A healthier, safer and more dignified walking experience.



Rollator Walker Buying Guide



Weight Capacity
(136 kg)



Total Weight
(9 kg)



Height range
71-98cm



Safety performance



Product Structure

- (1) Handle/brake
- (2) Tray
- (3) Seat
- (4) Lock lever

- (5) Storage basket
- (6) Item rack
- (7) Safety lock
- (8) Wheel

Product Advantage

Parking brake

Used as a recreational vehicle When you feel tired, you can lock, sit and rest. Pull the handle to lock the car and keep the car stable.

Brake deceleration

Allowing for better control and added safety. Feather-touch brakes offer easy manoeuvrability.

Folding structure, lightweight and portable

Saving space; it is very easy to load or store this Rollator into your vehicle, trains and airplanes.