

Rollator Walker

The rollator restores activity and independence for the elderly and mobility challenged. Is a light and versatile rollator which is approved for indoor and outdoor use, With a width of only 27", the Compact rollator is even easier to manoeuvre in tight spaces, allowing you to go where you want when out on the town or when travelling. A healthier, safer and more dignified walking experience.



Rollator Walker Buying Guide



Weight Capacity
(136 kg)



Total Weight
(9 kg)



Height range
80-96cm



Safety performance



Product Structure

- (1) Handle/brake
- (2) Backrest
- (3) Height adjustment knob
- (4) Cushion
- (5) Shopping basket
- (6) Wheel

Product Advantage

Widened and thickened seat

Suitable for different people, comfortable and breathable, The ergonomic shape and high density foam cushion with laminated memory foam on top provides unrivalled seating comfort.

Backrest Support

The Back is simply the most comfortable rigid rollator backrest available today. It features an ergonomically shaped soft padded cushion made from high density foam.