

# *Manual and Electric Nursing Bed*



Adjustable Hospital Beds for Home Care or assisted living environments are used to better meet the needs of loved ones who have reduced , or require special care due to a medical condition. Adjustable hospital beds offer a variety of positions to assist with care, eating, reading, washing hair, washing feet, defecationg, etting in and out of bed, and improve sleeping.

## *Electric Bed Buying Guide*



Weight Capacity



Functions and features



Material



Accessories



## *Product Structure*

- Headboard and footboard
- Guardrail
- Mattress
- Infusion stand
- Control Handle
- Motor
- Bedpan
- Casters (with lock)

## *Product Advantage*

- Flexibility of adjustment: Electric powered hospital beds are adjustable beds that offer recovering patients to manage the bed position .
- Comfort: If you pick a bed that is hassle-free and convenient to use, then a full electric bed will really be an exceptional choice for you.
- Mobility:The hospital beds is with premium castors and central locking braking system.

Bed accessories

### **1. Backrest up & down**

It is easy for patient to sit up down and do some daily life activities easily through electric remoter controller, lessening nursing daily routine for both patients and caregivers.

### **2. knee up & down**

With the electric remote controllers, patients can realize the leg up down to make the blood and pain circulation smoothly.

### **3. Washing head hair**

The shampoo basin is placed at the backrest of position, it is very useful and helpful for the paralyzed patients to wash head hair and feet.

### **4. Defecation function**

Provide convenient toilet functions for users who cannot take care of themselves.

### **5. Turn left and right**

This electric hospital bed can adjust the position of left and right via the remote control, it is convenient for the paralyzed patient. With this function, suffers can make the body blood flow smooth

### **6. Transfusion drip holes**

Standard accessories.