

Powerful *Electric* *Wheelchair*



New Lightweight foldable electric power wheelchair is preferred by everyone. It is Very Convenient as a Mobility Wheelchair, with Powerful Motors and High-Volume Battery, New Comfy Go Gives You Freedom. With the New Anti-Leaning Rear Design, It is Perfectly Stable and Safe.

Electric Wheelchair Buying Guide



500W Dual Motor
(250W * 2)



Lead-acid batteries
12km battery life



Removable Controller



Humanized design
Powerful and Foldable



Product Structure

- Joystick
- Back cushions
- Seat belt
- Handbrake
- Driving Mode Switch
- Leggings
- Footstool
- Wheel

Product Advantage

Power Wheelchair Features

- 1. Durable and Stylish
- 2. Travel Distance up to 12km.
- 3. EABS brake system,ramp packing more safe.
- 4. Folding structure, lightweight and portable,Saving space.
- 5. 360°smart controller,operate more sensitive.
- 6. Removable seat cushion,breathable and comfort.

Folding structure, lightweight and portable

- 1.Saving space; it is very easy to load or store this wheelchair into your vehicle, trains and airplanes.
- 2.Fully Assembled & Ready to Use: The wheelchair comes already assembled.



Product Comparison



| Model | Alphard Series 74501 | Alphard Series 74502 |
|------------------|----------------------------------|-----------------------------------|
| Total Weight | 24kg+3kg (battery packs) | 36kg+9kg (battery packs) |
| Bearing capacity | 120kg | 120kg |
| Charging Time | 6-8 hours | 6-8 hours |
| Max Speed | 6 km/h | 6 km/h |
| Range | 15km | 10-12km |
| Battery | Lithium Battery DC 24V 13 AH | Lead-acid batteries DC 24V 12AH |
| Wheels | 8"PU tyres * 12" Pneumatic tyres | 10"PU tyres * 16" Pneumatic tyres |