

# **Tagatustos ng commode** **wheel chair ng China**







# PORTABLE COMMODE WHEELCHAIR

Folding, Portable, with Wheels and Drop-Arm  
Wheeled Commode with Footrests



## PRODUCT STRUCTURE



## PRODUCT FEATURE



## PRODUCT ADVANTAGES



350lb capacity



Easily removable tool-free



Chrome-plated-coated steel tubing frame



Swivel casters 2 with locks

### WELDED CONSTRUCTION

Chrome-plated-coated steel tubing frame with welded construction. Tool-free assembly. also available in silver vein finish



### STANDARD PLASTIC TOILET SEAT

Injection-molded large plastic toilet seat, 12 quart bucket and cover come standard. bucket slides easily in and out of place



### EASY-TO-RELEASE ARM

Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.



### SWING-AWAY AND DETACHABLE FOOTRESTS

The chair also has humanized supportive shown with footrests. It is detachable and easy for you to clean or take away for easy standing up.



## PRODUCT PARAMETER

Parameter	Category
Series	Rigel Series
Model No.	10701
Width	18" between arms   20" Outside
Seat Dimensions	17.5" (W) x 17" (D)
Height	19.25"
Weight	29 lbs.
Carton Shipping Weight	34 lbs.
Weight Limit	350 lbs.