

Safe and easy solution for helps a person from room to room.

Feature

- Height adjustable handles to suit user requirements
- Height, width and angle adjustable PU leg supports, Padded lower leg supports can be adjusted tool free
- Low profile base facilitates user feet placement
- Sturdy steel frame provides stable support for user
- Handle has multiple gripping points and is angled to enable nose-over-toe standing movement
- Opening on central handle bar allows quick fitting of optional support belts
- 2 rear 8cm castors ensures effortless manoeuvring even with heavier users
- Easy wipe-clean surfaces

Product Name	Patient Turner
Net Weight	16.5 kg
Weight Capacity	150 kg
Height	6.3"-7.1"
Pad Dimensions	6.9" * 4.9"
Function	Perfect for room to room transfers
Application	Home care, nursing home or hospital, etc.

Detailed Description



Padded lower leg supports can be adjusted in width, angle and height for user comfort (tool free).



Low profile base facilitates user feet placement, Substantial base area provides stability, with rear bar also allowing carer to counterbalance.



Sturdy steel frame, opening on central handle bar allows quick fitting of optional support belts.



Weight capacity 150kg, It is helps carers or attendants to transfer a user from one seat to another.