



Lifting from floor



Toileting/Bathing transfer



Bed-Chair transfer

Patient lift is an assistive device, suitable for transfer people such as Wheelchair Users, Bariatric Patient, Bedridden patients, Disabled, Handicap, lower limb fracture patients, elderly people with poor ability to walk to a wheelchair.

Transfers may be to and from beds, chairs, floor to bed, lateral transfers, bathing, and toileting. The equipment is compact in design and versatile, suitable for moving all types of patients. Lightweight and easy to disassemble, convenient for fast transportation and installation. Can make patients and caregivers feel safe, comfortable and stable.

Features:

- To charge, simply plug in unit. Battery does not need to be removed to be charged.
- Lift will not operate while plugged in. Ideal for home-care use.
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely.
- Base legs adjust easily and lock securely into open position with padded shifter handle.
- Extra wide coated cradle with 360 rotation with 4-point hook up.
- Weight Capacity 600 lbs, (Can fit under bed).
- Optional Weight Scale.

Product Name	Bariatric Battery-Powered Lift
Net Weight	60kg
Weight Capacity	272kg
Height	32"-75.5"
Width of Base	25"-45.75"
Function	Handicapped or disabled transfer, walking assist, home care
Application	Home care, nursing home or hospital, etc.
Thrust	12000N

Detailed Description:



Easy operating

During the lifting process, the hanger can be rotated 360 to obtain a comfortable position.
Optional weight scale.



Emergency stop button

with Emergency stop button and emergency lowering system.
The pump can be switched to manual mode to lower the patient safely.



Adjustable base

Expand base of lifter to widest possible position to maximize stability and ease of use.



Castors with brake

Brakes in the rear wheels prevent the device from moving during use.