

This is a quiet and hard-working compressor nebulizer. It offers higher flow flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

**Detailed Description:**  
Ultra compact design makes it extremely portable.  
Powerful piston pump is efficient, dependable and quiet.  
Delivers consistent particle size of 4.5 microns and fast treatment times.



Suitable for A Variety of Situations




**PORTABLE  
COMPRESSOR NEBULIZER**

Penetrate into Bronchus & Lung 2-5µm



**MICRON NEBULIZATION  
FINE PARTICLES**

MMAD: 3.6 - 4.0 µm




**QUIET ENVIRONMENT  
SILENT ATOMIZATION**

Low noise design, easy for children to sleep



**USE BASIC STEPS**



More Details

