

## Battery Powered Patient Lift

This Lift offers exceptional value in an Electric Patient Lifter. It's an affordable powered lift with battery and charger. Plugs into any AC socket to charge, charger has a function light so you know when it is fully charged.



### Patient Lift Buying Guide



Weight capacity 204kg



Emergency lowering system



2 castors with brake brakes



Optional Weight Scale



## ***Product Structure***

- Emergency stop button
- Emergency lowering system
- Cradle
- Control Box
- Hand Control
- Push Handle
- Motor
- Spreader Bar
- Base
- Casters

P

## ***Product Advantage***

1. Easy to use Lift system has a wide cradle that accommodates a 2-point or 4-point sling.
2. Lift from floor, bed, commode, or wheelchair.
3. Heavy-gauge construction, with 6-Point Cradle.
4. Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely.
5. Fast, easy charge from AC outlet.
6. To charge, simply plug in unit. Battery does not need to be removed to be charged.
7. Lift will not operate while plugged in.
8. For home-care use.
9. Optional Weight Scale.