

Bariatric Battery-Powered Lift

Model 71922 Bariatric Electric Patient Lift is safe, stable and secure, giving users a feeling of confidence when being transferred. The patient lift has a 600-pound weight capacity, this bariatric battery-powered patient lift is equipped with a high-performance rechargeable battery and power to spare.



Patient Lift Buying Guide



Weight capacity 272kg



Emergency
Stop system



2 castors
with brake levers



Optional Weight Scale



Product Structure

- Cradle
- Mast & Boom
- Push Handle
- Motor
- Control Box
- Hand Control
- Spreader Bar
- Casters
- Base

Product Advantage

1. To charge, simply plug in unit Battery does not need to be removed to be charged
2. Lift will not operate while plugged in Ideal for home-care use
3. Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
4. Base legs adjust easily and lock securely into open position with padded shifter handle
5. Extra wide coated cradle with 360° rotation with 4-point hook up
6. Weight Capacity 600 lbs, (Can fit under bed)
7. Optional Weight Scale.