

# Sit to Stand Patient Mover 72130

## SPECIFICATIONS

Size(W*D*H)	22.6" * 24.8" * 46.5"
Base Plate Size	19.3" * 24.8"
Knee Pad Size	6.3"*7.1"
Leg Support Height	11.8-18.5"
Base Height	2"
Weight	36 lbs
Weight Limit	331 lbs




**331 lbs**  
Weight capacity

## DETAILED DESCRIPTION



A less expensive alternative to powered patient lifts, these standing assist devices are perfect for those capable of pulling themselves into a semi-standing posture and possessing enough leg strength to support themselves.

Maximum user weight of 150kg





Once the standing hoist is rolled into position, the patient simply grabs the center bar and hoists their body into a standing position. Opening on central handle bar allows quick fitting of optional support belts, and the handle has multiple gripping points.



Low ground clearance provides maximum stability during use.

The wheels feature easily engaged brakes and allow for movement over short distances; it's especially helpful during wheelchair transfers and toileting.



Adjustable knee pad, Contoured, height-adjustable  
knee pad offers high-level comfort and accommodates  
a large range of patient sizes.

---

## **SIMILAR PRODUCTS**

