

Procyon compressor nebulizer provides fast, efficient and reliable nebulized inhalation therapy, designed for home use. This is a quiet and hard-working compressor nebulizer suitable for nebulizing most aerosol medications prescribed for common respiratory diseases. It offers higher liter flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

Detailed Description: Makes it extremely portable.

It is compact, quiet and efficient.

Powerful, gentle, quiet & efficient nebulizer and quiet.

Delivers consistent particle size of 4.5 microns and fast treatment times.

Series / Model No.	Procyon series / 70000
Unit Dimensions	191mm x 102mm x 117mm
Weight	4.1kg (11.04kg)
Electrical Requirements	100-240VAC, 50/60Hz
Operating Noise Level	< 60 dBA
Maximum Pressure	< 15psi
Compressor flow air flow	> 6L/min
Particle Size	2.0 - 5.0µm
Nebulizer Flow	> 2.0L/min
Residual Liquid	10.5ml




**PORTABLE
COMPRESSOR NEBULIZER**

Penetrate into Bronchus & Lung 2-5µm



**MICRON NEBULIZATION
FINE PARTICLES**

MMAD: 3.6 - 4.0 µm




**QUIET ENVIRONMENT
SILENT ATOMIZATION**

Low noise design, easy for children to sleep




USE BASIC STEPS



Suitable for A Variety of Situations



More Details



Micron level administration—2µm–5µm
Nebulized particles directly reach the upper respiratory tract, trachea, bronchi and lower respiratory tract, making it easier to absorb the medication.



High atomization efficiency
Benefits atomization for convenience operation.



Low noise design
Low power, low vibration, baby atomization does not exist.