

Procyon compressor nebulizer provides fast, efficient and reliable nebulized inhalation therapy, designed for home use. This is a quiet and hard-working compressor nebulizer suitable for nebulizing most aerosol medications prescribed for common respiratory diseases. It offers higher liter flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

Performance and Design Features:

- Reliable pressure and flow
- Reusable nebulizer can be included
- Small footprint
- Small, compact design
- Convenient carrying handle
- Metal mechanical parts

Weight: 2.18kg



Series / Model No.	Procyon series / 70802
Unit Dimensions	191mm*173mm*102mm
Weight	2.18kg
Electrical Requirements	100-240VAC, 50-60HZ
Operating Noise Level	≤ 60 dBA
Maximum Pressure	≥ 35psi
Compressor free air flow	≥ 6L/min
Particle Size	2µm - 5µm
Nebulizer Rate	≥ 0.2ml/min
Residual Liquid	70.5ml




PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



Adults & Children

MICRON NEBULIZATION FINE PARTICLES



<p>3.6-4.0µm Median particle size MMAD</p>	<p>75% Percentage of effective particles</p>
---	---

QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



< 60dB
Noise Level!

Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

