

This is a quiet and hard-working compressor nebulizer. It offers higher liter flows, shorter treatment times and lower sound levels than many other economy compressor nebulizers.

Detailed Description:


Ultra compact design makes it extremely portable.

Powerful piston pump is efficient dependable and quiet.


Delivers consistent particle size of ≤5 microns and fast treatment times.



Suitable for A Variety of Situations



More Details




Micron level nebulization—0.5 μ m~5 μ m

Nebulized particles directly reach the upper respiratory tract, trachea, alveoli and lower respiratory tract, making it easier to absorb the medication.



High atomization efficiency

Smooth atomization for continuous operation.



Low noise design

Low power, low vibration, baby atomization does not resist.

