

## Feature

- Opening on central handle bar allows quick fitting of optional support belts
- 4 front and rear castors combined with 2 x 13cm (5") central wheels ensure effortless manoeuvring even with heavier users
- Foot operated brakes on centre wheels provide stability during standing/sitting movement
- Easy wipe clean surfaces
- Easily dismantled/reassembled without tools for transportation and storage, it has a Maximum User Weight of 150kg
- Substantial base area provides stability

<b>Product Name</b>	Sit to Stand Transfer Aid
<b>Dimensions</b>	22.6" * 24.8" * 46.5"
<b>Weight Capacity</b>	150 kg
<b>Height</b>	11.8" - 18.5"
<b>Function</b>	Perfect for room to room transfers
<b>Application</b>	Home care, nursing home or hospital, etc.

## Detailed Description



Sturdy steel frame, opening on central handle bar allows quick fitting of optional support belts.



Padded lower leg support can be adjusted in width, angle and height to suit user comfort.



4 front and rear castors combined with 2 x 13cm (5") central wheels ensure effortless manoeuvring even with heavier users. it has a Maximum User Weight of 150kg.