

# Hydraulic Patient Lift

The Deneb series lifts are very easy to maneuver , and offers a wide selection slings for a variety of lifting situations. Transfers may be to and from beds, chairs, floor to bed, lateral transfers, bathing, and toileting .The equipment is compact in design and versatile, suitable for moving all types of patients. Lightweight and easy to disassemble,convenient for fast transportation and installation. Can make patients and caregivers feel safe, comfortable and stable.



## Patient Lift Buying Guide



Weight Capacity  
(450lbs)



Lift Sizing Choose



Six-point swivel bar



Hydraulic pump



### Product Structure

- Hanger
- Sling
- Base
- Rear wheel
- Base adjustment lever
- Hydraulic pump
- Boom
- Push hand

### Product Advantage

#### Patient lift ([Hoyer lift](#)) Features

- The Lift features a 20"-64" range and allows individuals to be picked up from lying position on floor
- Pump handle: able to rotate from side to side for caregiver convenience. How many times can be used.
- Tips: It is highly recommended that patients are assisted by two people while lifting.
- Easy Assembling, It Only takes you 2 minute to assembling the lifts.

#### [Lift Sizing Choose](#)



### Sling Features

Easy to use: Sling aids in transferring users from a bed to a wheelchair, recliner or shower chair. The large opening in the design also accommodates toileting.

Maintenance(Easy to clean): Washable polyester fabric can be used in wet or dry environments.

Multiple Use: Ideal for transferring a patient from the bed to a chair, wheelchair or toilet, chair to chair or from the floor to a bed or chair using your manual or electric lift.