

Patient Turner

Patient Turner, Safe and easy solution for transferring from one seated position to another, Easily dismantled/reassembled without tools for transportation and storage.

The Ablestand Patient Turner helps

1. Attendant to assist a patient
2. User to transfer to another seat, bed, wheelchair, toilet or commode by standing and turning



Patient Turner Buying Guide



Weight capacity 150kg



Sturdy steel frame



Height adjustable handles



8cm castors*2



Product Structure

- Adjustable handles
- Sturdy steel frame
- Knob
- Padded lower leg supports
- 3 "castors
- Low profile base

Product Advantage

- Height, width and angle adjustable PU leg supports,Padded lower leg supports can be adjusted tool free.
- Sturdy steel frame provides stable support for user
- Handle has multiple gripping points and is angled to enable nose-over-toe standing movement
- Easy wipe-clean surfaces
- Weight capacity 150kg
- User will need to pull on handle to stand and be weight bearing