

O nebulizador compressor Procyon fornece terapia de inalação nebulizada rápida, eficiente e confiável, projetada para uso doméstico. Este é um nebulizador compressor silencioso e resistente, adequado para nebulizar a maioria dos medicamentos em aerosol prescritos para doenças respiratórias comuns. Ele oferece fluxos de litros mais altos, tempos de tratamento mais curtos e níveis de ruído mais baixos do que muitos outros nebulizadores de compressores econômicos.

**Recursos de desempenho e design:**

- Pressão e fluxo confiáveis
- Nebulizador reutilizável pode ser incluído
- Pequena pegada
- Design pequeno e compacto
- Alça de transporte conveniente
- Peças mecânicas metálicas

Weight: 2.18kg



Série/modelo não.	Série Procyon / 70802
Dimensões da Unidade	191mm*173mm*102mm
Peso	2.18kg
Requisitos elétricos	100-240 VCA, 50-60 Hz
Nível de ruído operacional	≈ 60dBx
Pressão Máxima	≈ 35 psi
Fluxo de ar livre do compressor	≈ 6L/min
Tamanho da partícula	2µm - 5µm
Taxa do nebulizador	≈ 0.2 ml/min
Líquido Residual	<0.2ml




## PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



## MICRON NEBULIZATION FINE PARTICLES



<b>3.6-4.0µm</b> Median particle size MMAD	<b>75%</b> Percentage of effective particles
--	--

## QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

