

Fornecedor de cadeira de rodas de cômoda da china

PORTABLE COMMODE WHEELCHAIR

Folding, Portable, with Wheels and Drop-Arm
Wheeled Commode with Footrests



PRODUCT STRUCTURE



PRODUCT FEATURE



Tool-free removable(back)



Rear wheel
Front wheel
Rust-resistant, swivel casters
(2 with locks)

WELDED CONSTRUCTION

Chrome-plated-coated steel tubing
frame with welded construction.
Tool-free assembly.
also available in silver vein finish



STANDARD PLASTIC TOILET SEAT

Injection-molded large plastic
toilet seat, 12 quart bucket and
cover come standard.

bucket slides easily in and out of place



EASY-TO-RELEASE ARM

Easy-to-release arm mechanism allows
for safe lateral patient transfers to and
from commode.



SWING-AWAY AND DETACHABLE FOOTRESTS

The chair also has humanized supportive
shown with footrests.It is detachable
and easy for you to clean or take away
for easy standing up.



PRODUCT ADVANTAGES



350lb capacity



Easily removable
tool-free



Chrome-plated-coated
steel tubing frame



Swivel casters
2 with locks

PRODUCT PARAMETER

Parameter	Category
Series	Rigel Series
Model No.	10701
Width	18" between arms 20" Outside
Seat Dimensions	17.5" (W) x 17" (D)
Height	19.25"
Weight	29 lbs.
Carton Shipping Weight	34 lbs.
Weight Limit	350 lbs.