

# **China fabricante de** **cadeiras de rodas elétricas**







# TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



## PRODUCT ADVANTAGES

EABS Brake system	500W dual motor	High battery life	360° joystick
Hand/electric mode	Easy to fold	Aluminum alloy frame	Anti-backover design

## PRODUCT PARAMETER

Parameter	Category
Series/Model No.	Alphard Series/74501
Frame material	Aluminum alloy
Unfolded Dimension	1000*610*930mm
Folded Dimension	660*380*770mm
Motor	DC 24V、250W
Battery	Lithium battery DC 24V 13 AH
Weight	24KG
Loading Capacity	120KG
Front wheel specifications	8"*2" (200*50mm) PU tire
Rear wheel specifications	12"*2" (305*60mm) pneumatic tire
Max Speed	6KM/H
Driving Range	15KM

## PRODUCT STRUCTURE



## PRODUCT DETAILS



### FOLDING DESIGN

Lightweight folding portable design  
Can be placed in the trunk, Similar to suitcase size.



### LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



### THE ARMREST CAN BE RAISED

Handrails can swing up and down  
Easy to hold the patient up and down wheelchair



### HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel

