

Detailed Description:
Ultra compact design makes it extremely portable.
Powerful jetting pump is efficient, dependable and quiet.
Delivers consistent particle size of 45 microns and fast treatment times.



Suitable for A Variety of Situations




PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



MICRON NEBULIZATION FINE PARTICLES

MMAD: 3.6 - 4.0 µm




QUIET ENVIRONMENT SILENT ATOMIZATION


Low noise design, easy for children to sleep



USE BASIC STEPS



More Details




10-15µm Upper Airway

5-10µm Air ducts

2-5µm Lower Airway


Micro level subinhalation——2µm-5µm
Nebulized particles directly reach the upper respiratory tract, trachea, alveoli and lower respiratory tract, making it easier to absorb the medication.



Compressor free air flow: 6L/min

Nebuliser Rate: ≥0.2ml/min

High atomization efficiency
Smooth atomization for continuous operation.



<60 dBA

Low noise design
Proven, low vibration, baby atomization does not resist.