

□□□□□□□□ □□□□□□□□□□□□□□□□

□□□□□□□□□□□□□□□□ **commode**

PORTABLE COMMODE WHEELCHAIR

Folding, Portable, with Wheels and Drop-Arm
Wheeled Commode with Footrests



PRODUCT STRUCTURE



PRODUCT FEATURE



PRODUCT ADVANTAGES



350lb capacity



Easily removable tool-free



Chrome-plated-coated steel tubing frame



Swivel casters 2 with locks

WELDED CONSTRUCTION

Chrome-plated-coated steel tubing frame with welded construction. Tool-free assembly. also available in silver vein finish



STANDARD PLASTIC TOILET SEAT

Injection-molded large plastic toilet seat, 12 quart bucket and cover come standard. bucket slides easily in and out of place



EASY-TO-RELEASE ARM

Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.



SWING-AWAY AND DETACHABLE FOOTRESTS

The chair also has humanized supportive shown with footrests. It is detachable and easy for you to clean or take away for easy standing up.



PRODUCT PARAMETER

Parameter	Category
Series	Rigel Series
Model No.	10701
Width	18" between arms 20" Outside
Seat Dimensions	17.5" (W) x 17" (D)
Height	19.25"
Weight	29 lbs.
Carton Shipping Weight	34 lbs.
Weight Limit	350 lbs.